

高.低溫攪拌恆溫槽

Constant Temp. Bath with Magnetic Stirrers

PSL-1400 · 1810 · 2000

高溫合成用

2 L
室溫 ~ 150°C



PS-1000

300 mL
-40 ~ 0°C



PSL-1400

低溫合成用

300 mL
-80 ~ 0°C



PSL-1810

2 L
-80 ~ 20°C




PSL-2000




PSL-1400 操作面板




PSL-1810 操作面板




PSL-2000 操作面板



PSL-1000




PSL-1400 · 1810



PSL-2000



安全固定夾



可調式上蓋 (PSL-1800付)



結晶研究
(PSL-2000)

型 號	PS-1000	PSL-1400	PSL-1810	PSL-2000
溫度調節範圍	室溫 +10 ~ 150°C (油) 室溫 +5 ~ 80°C (水)	-40 ~ 0°C	-80 ~ 0°C	-80 ~ 20°C
溫度調節精度	± 0.5°C (水)、± 1.5°C (油)	± 0.5°C 以內	± 1.0°C 以內	± 1.0°C 以內
轉速範圍	100 ~ 800 rpm			
加熱功率	700 W	250 W	450 W	750 W
冷凍機	-	350 W	450 W × 2	500 W × 2
內槽尺寸：容量	φ 250 × 130 D mm, 6 L	φ 140 × 135 D mm, 2.1 L		φ 250 × 180 D mm, 8 L
可用容器量大容量 / 直徑	2 L, φ 165 mm	300 mL, φ 95 mm		2 L, φ 135 mm

* 另有可放 9 支 φ 10 ~ 18mm 試管用 holder 可供選擇