

# 切割式粉碎機

Cutting MILL/Cross Beater Mill /Rotor Beater Mill



SM2000



SM100



SK100



SR200/SR300

## 規格說明

| 型號      | SM100                            | SM2000     | SK100               | SR200               | SR300      |
|---------|----------------------------------|------------|---------------------|---------------------|------------|
| 適用材質    | Soft、medium-hard、elastic、fibrous |            | medium-hard、brittle | soft to medium-hard |            |
| 入料尺寸    | Max.60×80 mm                     |            | 15 mm               | 15 mm               | 15 mm      |
| 出料尺寸    | 0.25~20 mm (視網篩而定)               |            | <100 μm             | <80 μm              | <50 μm     |
| 入料方式    | Continuous                       | Continuous | Continuous          | Continuous          | Continuous |
| 最大樣品量   | 5 or 30 L                        | 5 or 30 L  | 5 or 30 L           | 5 or 30 L           | 5 or 30 L  |
| Sieve尺寸 | 0.25~20 mm                       | 0.25~20 mm | 0.12~10 mm          | 0.08~10 mm          | 0.08~10 mm |
| Rotor材質 | 不鏽鋼及 Heavy-metal Free steel 可選   |            | 不鏽鋼、鉻鋼              | 鑄鐵                  | 不鏽鋼        |

# 刀片式/旋轉式粉碎機

Knife / Rotor Mill



GM200



ZM200

## 規格說明

| 型號    | GM200                       | 型號       | ZM200                                |
|-------|-----------------------------|----------|--------------------------------------|
| 適用材質  | Soft、medium-Hard、Elastic 材料 | 適用材質     | Soft、Medium-Hard、Elastic、Fibrous 材料  |
| 入料尺寸  | 10~40 mm                    | 入料尺寸     | 10 mm                                |
| 出料尺寸  | <300 μm                     | 出料尺寸     | <40 μm                               |
| 入料方式  | Batch                       | 入料方式     | Continuous                           |
| 研磨杯材質 | PP塑膠、PC塑膠、強化玻璃、不鏽鋼 可選       | Rotor 材質 | 不鏽鋼、不鏽鋼包覆碳化鎢、鈦金屬 可選                  |
| 研磨杯大小 | 500ml、1000 ml 可選            | Sieve 尺寸 | 0.08 - 6 mm                          |
| 最大樣品量 | 300ml - 700 ml              | Sieve 規格 | 6 tooth、8 tooth、12 tooth、24 tooth 可選 |
| 轉速設定  | 2000 - 1000 rpm             | 最大樣品量    | 900 ~ 5000 ml                        |
| 時間設定  | 1sec - 3 min                | 轉速設定     | 6000 ~ 18000 rpm                     |