


8爪及12爪微量吸管

Eppendorf Research® multichannel pipettes



提供樣品儲存試劑槽

不同的吸取範圍選擇

Description

Long series in columns, rows and plates—that is the area in which the multichannel pipette is becoming increasingly popular. And it is especially for this environment that we developed the **Research multichannel** pipette, which can be operated easily with one hand. This high-performance 8-or 12-channel pipette offers an effective, precise and rapid way to process long series for an entire spectrum of applications.



Each Eppendorf Research multichannel pipette contains a service kit and an epTIPS Box with the appropriate volume tip.

Ordering information : Eppendorf Research multichannel pipettes

Volume range	Description	Volume	Volume increment	Inaccuracy*	Imprecision*	Order no.
gray operating button, for use with 20 µl pipette tips						
0.5 - 10 µl	8 - channel	1 µl	0.01 µl	± 8.0 %	≦ 5.0 %	3114 000.115
	12 - channel	5 µl 10 µl		± 4.0 % ± 2.0 %	≦ 2.0 % ≦ 1.0 %	
yellow operating button, for use with 200 µl pipette tips						
10 - 100 µl	8 - channel	10 µl	0.1 µl	± 3.0 %	≦ 2.0 %	3114 000.131
	12 - channel	50 µl 100 µl		± 1.0 % ± 0.8 %	≦ 0.8 % ≦ 0.3 %	
yellow operating button, for use with 300 µl pipette tips						
30 - 300 µl	8 - channel	30 µl	0.2 µl	± 3.0 %	≦ 1.0 %	3114 000.158
	12 - channel	150 µl 300 µl		± 1.0 % ± 0.6 %	≦ 0.5 % ≦ 0.3 %	

Ordering information : Accessories

Description	Order no.
Lower parts for:	
0.5 - 10 µl 8 - channel	3114 608.008
12 - channel	3114 609.004
10 - 100 µl 8 - channel	3114 610.002
12 - channel	3114 611.009
30 - 300 µl 8 - channel	3114 612.005
12 - channel	3114 613.001
Reagent Reservoir , autoclavable reservoir for reagent takeup with multichannel pipette, filling volume 50 - 60 ml, 1 set = 10 pcs.	0030 058.607
Pipette Carousel Stand , for 6 pipettes	3115 000.003